

Anti-BICD2 antibody (440-620) (STJ112330)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 94kDa Observed Mw: 118kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23299
Gene Symbol	BICD2
UniProt ID	BICD2_HUMAN
Immunogen Region	440-620
Immunogen Sequence	HVAVAEAGELREQLKALRST HEAREAQAEEKGRYEAEGQ ALTEKVSLLLEKASRQDRELL ARLEKELKKVSDVAGETQGS LSVAQDELVTTFSEELANLYH HVCMCNNETPNRVMLDYRE GQGGAGRTSPGGRTSPEARG RRSPIILLPKGLLAPEAGRAD GGTGDSSSPSGSSLPSPLSD P
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 440-620 of human BICD2 (NP_001003800.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using BICD2 antibody (STJ112330) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.