

Anti-MIEF1 antibody (40-280) (STJ112432)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 51kDa Observed Mw: 51kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54471
Gene Symbol	MIEF1
UniProt ID	MID51_HUMAN
Immunogen Region	40-280
Immunogen Sequence	GIATLAVKRMVDRAISAPTS PTRLSHSGKRSWEEPNNWMS PRLNDRDMKTGLSRSLQTLPTDSSTFDTDTCPPRPKPVA RKGQVDLKKSLRMSLQEKLTYYRNRAAIPAGEQARAKQAAVDICAE LRSFLRAKL PDM PLRDMYLSGSLYDDLQVVTA DHIQLIVPLVLEQNLWSCIPGEDTIMNVPGFFLVRRNPE YFPRGSSYWDRCVVGGYLSP KTVADTFEKVVAGSINWPA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 40-280 of human MIEF1 (NP_061881.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates, using MIEF1 antibody (STJ112432) at 1:400 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.