

## Anti-NOXO1 antibody (121-370) (STJ112503)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 $\mu$ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 41kDa Observed Mw: 41kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">124056</a>
Gene Symbol	<a href="#">NOXO1</a>
UniProt ID	<a href="#">NOXO1_HUMAN</a>
Immunogen Region	121-370
Immunogen Sequence	QPLDLEPALPPGSRVILPTP EEQPLSRAAGRLSIHSLEAQ SLRCLQPFCTQDTRDRPFQA QAQESLDVLLRHPSGWWLVE NEDRQTAWFPAPYLEEAAPG QGREGGPSLGSSGPQFCASR AYESSRADELSVPAGARVRV LETSDRGWWLCRYGDRAGLL PAVLLRPEGLGALLSGTGFR GGDDPAGEARGFPEPSQATA PPPTVPTRPSPGAIQSRCT VTRRALERRPRRQGRPRGC
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 121-370 of human NOXO1 (NP_653204.1).

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of extracts of various cell lines, using NOXO1 antibody (STJ112503) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25  $\mu$ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.