

Anti-CSF2RB antibody (200-443) (STJ112525)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 97kDa Observed Mw: 95-120kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1439
Gene Symbol	CSF2RB
UniProt ID	IL3RB_HUMAN
Immunogen Region	200-443
Immunogen Sequence	EHLMPSSTYVARVTRLAPG SRLSGRPSKWSPEVCWDSQP GDEAQPQNLECFDGAAVLS CSWEVRKEVASSVSFGLFYK PSPDAGEEECSVLRGLGS LHTRHHCQIPVDPATHGQY IVSVQPRRAEKHIKSSVNIQ MAPPSLNVTKDGDYSYLRWE TMKMRYEHIDHTFEIQYRKD TATWKDSTETLQNAHSMAL PALEPSTRYWARVRVRTSRT GYNGIWSEWSEARSWDTES
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 200-443 of human CSF2RB (NP_000386.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates, using CSF2RB Rabbit polyclonal antibody (STJ112525) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.