

Anti-COPG1 antibody (555-874) (STJ112569)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 98kDa Observed Mw: 105kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22820
Gene Symbol	COPG1
UniProt ID	COPG1_HUMAN
Immunogen Region	555-874
Immunogen Sequence	IPGLERALQQYTLEPSEKPF DLKSVPLATAPMAEQRTTEST PITAVKQPEKVAATRQEIFQ EQLAAVPEFRGLGPLFKSSP EPVALTESETEYVIRCTKHT FTNHMVVFQFDCTNTLNDQTL ENVTVQMEPTAEAYELCYVP ARSLPYNQPGTCYTLVALPK EDPTAVACTFSCMMKFTVKD CDPTTGETDDEGYEDEYVLE DLEVTVDHIQKVMKLNFEA AWDEVGDFEKEETFTLST
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 555-874 of human COPG1 (NP_057212.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates using COPG1 Rabbit polyclonal antibody (STJ112569) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.