

Anti-GUCY2C antibody (24-330) (STJ112579)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 123kDa Observed Mw: 123kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2984
Gene Symbol	GUCY2C
UniProt ID	GUC2C_HUMAN
Immunogen Region	24-330
Immunogen Sequence	SQVSNCHNGSYEISVLMMG NSAFAPLKNLEDAVNEGLE IVRGRLQNAAGLNVTVNATFM YSDGLIHNSGDCRSSTCEGL DLLRKISNAQRMGCVLIGPS CTYSTFQMYLDTLSYPMIS AGSFGLSCDYKETLTRLMSP ARKLMYFLVNFWKTNLDPFK TYSWSTSYYKNGTETEDCF WYLNALASVSYSFSELHGFK VVLRQDKEFQDILMDHNRKS NVIIMCGGPEFLYKLGDR
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 24-330 of human GUCY2C (NP_004954.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of various cell lines, using GUCY2C antibody (STJ112579) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.