

Anti-WDR11 antibody (975-1224) (STJ112603)

GENERAL INFORMATION

	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 137kDa Observed Mw: 130kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55717
Gene Symbol	WDR11
UniProt ID	WDR11_HUMAN
Immunogen Region	975-1224
Immunogen Sequence	YFQKFLERVNLQEVKRSTY DHTRKCTDQLLLLGQTDRAV QLLLETSADNQHHYCDLSKA CLVTTVTSSGSPSQSTIKLVA TNMIANGKLAEGVQLLCLID KAADACRYLQTYGEWNRAAW LAKVRLNPEECADVLRRWVD HLCSPQVNNQKSKALLVLLSL GCFFSVAETLHSMRYFDRAA LFVEACLKYGAFEVTEDETEK LITAIYADYARSLKNLGFQK GAVLFASKAGAAGKDLLNE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 975-1224 of human WDR11 (NP_060587.8).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of various cell lines, using WDR11 antibody (STJ112603) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.