

Anti-LIMA1 antibody (560-759) (STJ113282)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:100-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 85kDa Observed Mw: 100kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51474
Gene Symbol	LIMA1
UniProt ID	LIMA1_HUMAN
Immunogen Region	560-759
Immunogen Sequence	DEISKPEVPEDVDLDLKKLR RSSSLKERSRPFTVAASFQS TSVKSPKTVSPPIRKGWSMS EQSEESVGGGRVAERKQVENA KASKKNGNVGKTTWQNKESK GETGKRSKEGHSLEMENENL VENGADSDDEDDNSFLKQQSP QEPKSLNWSSFVDNTFAEEF TTQNQKSQDVELWEGEVVKE LSVEEQIKRNRYYDEDEDEE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 560-759 of human LIMA1 (NP_057441.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of lysates from A-549 cells, using LIMA1 Rabbit polyclonal antibody (STJ113282) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.