

## Anti-SPAG6 antibody (180-290) (STJ113564)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human              |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                       |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                            |
| Concentration       | Lot specific                                                                                                                                          |
| Conjugation         | Unconjugated                                                                                                                                          |
| Purification        | Affinity purification                                                                                                                                 |
| Dilution Range      | WB:1:500-1:1000<br>ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                                      |
| Isotype             | IgG                                                                                                                                                   |
| Molecular Weight    | Protein Mw: 55kDa<br>Observed Mw: 55kDa                                                                                                               |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                        |

### TARGET INFORMATION

|                    |                                                                                                                         |
|--------------------|-------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">9576</a>                                                                                                    |
| Gene Symbol        | <a href="#">SPAG6</a>                                                                                                   |
| UniProt ID         | <a href="#">SPAG6_HUMAN</a>                                                                                             |
| Immunogen Region   | 180-290                                                                                                                 |
| Immunogen Sequence | EIALKRIAASALSDIAKHSP ELAQTVVDAGAVAHQAQML NPDALKLKHQILSALSQVSKH<br>SVDLAEMVVEAEIFPVVLTG LKDKDEYVKKNASTLIREIA KHTPELSQLVV |
| Specificity        | Recombinant fusion protein containing a sequence corresponding to amino acids 180-290 of human SPAG6 (NP_758442.1).     |

### ADDITIONAL INFORMATION

**Note** STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using SPAG6 Rabbit polyclonal antibody (STJ113564).  
 Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution.  
 Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.