

Anti-SIX1 antibody (130-284) (STJ113870)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 32kDa Observed Mw: 36kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6495
Gene Symbol	SIX1
UniProt ID	SIX1_HUMAN
Immunogen Region	130-284
Immunogen Sequence	CFKEKSRGVLREWYAHNPYP SPREKRELAETGLTTTQVS NWFKNRRQRDRAAEAKEREN TENNSSSNKQNLSPLEGG KPLMSSSEEEFSPQSPDQN SVLLLLQGNMGHARSSNYSLP GLTASQPSHGLQTHQHQLQD S
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 130-270 of human SIX1 (NP_005973.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of C2C12, using SIX1 antibody (STJ113870) at 1:800 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.