

## Anti-RNF4 antibody (1-190) (STJ113938)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Mouse/Rat          |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                           |
| Concentration       | Lot specific                                                                                                                                         |
| Conjugation         | Unconjugated                                                                                                                                         |
| Purification        | Affinity purification                                                                                                                                |
| Dilution Range      | WB:1:500-1:1000<br>ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                                     |
| Isotype             | IgG                                                                                                                                                  |
| Molecular Weight    | Protein Mw: 21kDa<br>Observed Mw: 21kDa                                                                                                              |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                       |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                               |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">6047</a>                                                                                                                                                                                          |
| Gene Symbol        | <a href="#">RNF4</a>                                                                                                                                                                                          |
| UniProt ID         | <a href="#">RNF4_HUMAN</a>                                                                                                                                                                                    |
| Immunogen Region   | 1-190                                                                                                                                                                                                         |
| Immunogen Sequence | MSTRKRRGGAINSRQAQKRT REATSTPEISLEAEPIELVE TAGDEIVDLTCESELPVVVD<br>LTHNDSVVVDERRRRPRRNA RRLPQDHADSCVVSSDDEEL SRDRDVYVTTHTPRNARDEG<br>ATGLRPSGTVSCPICMDGYS EIVQNGRLIVSTECGHVFCS QCLRDSLKNANTCPTCRKKI NHKRYHPIYI |
| Specificity        | Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human RNF4 (NP_001171938.1).                                                                                           |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Western blot analysis of various lysates using RNF4 Rabbit polyclonal antibody (STJ113938) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.