

Anti-SEMA7A antibody (45-280) (STJ113940)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 75kDa Observed Mw: 75kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8482
Gene Symbol	SEMA7A
UniProt ID	SEM7A_HUMAN
Immunogen Region	45-280
Immunogen Sequence	QGHLRSGPRIFAVWKGHVQ DRVDFGQTEPHTVLFHEPGS SSVWVGGRGKVYLFDFPEGK NASVRTVNIGSTKGSCLDKR DCENYITLLERRSEGLLACG TNARHPSCWNLVNGTVVPLG EMRGYAPFSPDENSLVLFEG DEVYSTIRKQEYNGKIPRFR RIRGESELYTSDTVMQNPQF IKATIVHQDQAYDDKIYYFF REDNPDKNPEAPLNVSRVAQ LCRGDQGGESSLSVSK
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 45-280 of human SEMA7A (NP_003603.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates using SEMA7A Rabbit polyclonal antibody (STJ113940) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.