

Anti-EXOC4 antibody (1-280) (STJ114252)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 110kDa Observed Mw: 122kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	60412
Gene Symbol	EXOC4
UniProt ID	EXOC4_HUMAN
Immunogen Region	1-280
Immunogen Sequence	MAAEAAGGKYRSTVSKSKDP SGLLISVIRTLSTSDDEVDR ENEKGRLEEAYEKCDRLDE LIVQHYTELTTAIRTYQSIT ERITNSRNKIKQVKENLLSC KMLLHCKRDELRLKLVWIEGIE HKHVLNLLDEIENIKQVPQK LEQCMASKHYLSATDMLVSA VESLEGPLLQVEGLSDLRLE LHSKKMNLLHLVIDELHRHL YIKSTSRVVRNKEKGISS LVKDASVPLIDVTNLPTPR
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human EXOC4 (NP_068579.3).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates using EXOC4 Rabbit polyclonal antibody (STJ114252) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.