

Anti-SRSF9 antibody (1-221) (STJ114412)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 26kDa Observed Mw: 24kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	8683
Gene Symbol	SRSF9
UniProt ID	SRSF9_HUMAN
Immunogen Region	1-221
Immunogen Sequence	MSGWADERGGEGDGRIYVGN LPTDVREKDLLEDFYKYGRI REIELKNRHGLVPPAFVRFEDPRDAEDIAIYGRNGYDYGQC RLRVEFPRTYGGRRGWPRGG RNPPTRRSDFRVLVSGLPSSGSWQDLKDHMREAGDVCYA DVQKDGVMVEYLRKEDMEY ALRKLDDTKFRSHEGETSYIRVYPERSTSYGYSRSRSGSR GRDSPYQSRGSPHYFSPFRP Y
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-221 of human SFRS9 (NP_003760.1).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Immunofluorescence analysis of U2OS cells using SRSF9 antibody (STJ114412) at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of extracts of various cell lines, using SRSF9 antibody (STJ114412) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.