

## Anti-MAN1C1 antibody (44-240) (STJ114703)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/ELISA           |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse        |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                      |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                           |
| Concentration       | Lot specific                                                                                                                                         |
| Conjugation         | Unconjugated                                                                                                                                         |
| Purification        | Affinity purification                                                                                                                                |
| Dilution Range      | WB:1:500-1:2000<br>ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                                     |
| Isotype             | IgG                                                                                                                                                  |
| Molecular Weight    | Protein Mw: 71kDa<br>Observed Mw: 68kDa                                                                                                              |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                       |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                    |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">57134</a>                                                                                                                                                                                              |
| Gene Symbol        | <a href="#">MAN1C1</a>                                                                                                                                                                                             |
| UniProt ID         | <a href="#">MA1C1_HUMAN</a>                                                                                                                                                                                        |
| Immunogen Region   | 44-240                                                                                                                                                                                                             |
| Immunogen Sequence | PHSSRLKRLFLAPRTQQPGL EVVAEIAGHAPAREQEPPPN PAPAAPAGGEDDPSSWASPR<br>RRKGGLRRTRPTGPREETA ARGNSIPASRPGEVPRF DFNAFRSLRHPVLGTRADE<br>SQEPQSQVRAQREKIKEMMQ FAWQSYKRYAMGKNELRPLT KDGYESNMFGGLSGATVIDS<br>LDTLYLMELKEEFQEAQ |
| Specificity        | Recombinant fusion protein containing a sequence corresponding to amino acids 44-240 of human MAN1C1 (NP_065112.1).                                                                                                |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



Western blot analysis of various lysates using MAN1C1 Rabbit polyclonal antibody (STJ114703) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.