

Anti-PARP8 antibody (1-300) (STJ114760)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 96kDa Observed Mw: 96kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	79668
Gene Symbol	PARP8
UniProt ID	PARP8_HUMAN
Immunogen Region	1-300
Immunogen Sequence	MGMCSRQERIQKDIDVVIQK SRAEKDCLFADFRYSDSTFT FTYVGGPRSVSYSVHVSEDY PDNTYVSSSENEDVLVTTE PIPVIFHRIATELRKTNIDN CCLSIKSKLOKENGEEESRQN STVEEDSEGDNDSEEFYYGG QVNYDGELHKHPQLEADLSA VREIYGPHAVSLREYGAIDD VDIDLHIDVSFLDEEIAVAW EVIRTEPIIVRLHCSLTQYL NGPVPTVDVVFQISTKERFG
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human PARP8 (NP_001171527.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunohistochemistry analysis of PARP8 in paraffin-embedded Human placenta using PARP8 Rabbit polyclonal antibody (STJ114760) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using PARP8 Rabbit polyclonal antibody (STJ114760) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.