

Anti-PDIA2 antibody (326-525) (STJ114937)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 58kDa Observed Mw: 58kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	64714
Gene Symbol	PDIA2
UniProt ID	PDIA2_HUMAN
Immunogen Region	326-525
Immunogen Sequence	FGLKAEAPTLLRLVNLETTK KYAPVDGGPVTAAASITAFCH AVLNGQVKPYLLSQEIPPDW DQRPVKTLVGKNFEQVAFDE TKNVFKFYAPWCTHCKEMA PAWEALAEKYQDHEDIIIAE LDATANELDFAVHGFPTLK YFPAGPGRKVIEYKSTRDLE TFSKFLDNGGVLPTTEEPPEE PAAPFPEPPANSTMGSKKEEL
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 326-525 of human PDIA2 (NP_006840.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of NIH/3T3 cells using PDIA2 Rabbit polyclonal antibody (STJ114937) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using PDIA2 Rabbit polyclonal antibody (STJ114937) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.