


Anti-IFI30 antibody (1-250) (STJ114964)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
	Isotype	IgG
	Molecular Weight	Protein Mw: 28kDa Observed Mw: 28kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	10437
	Gene Symbol	IFI30
	UniProt ID	GILT_HUMAN
	Immunogen Region	1-250
	Immunogen Sequence	MTLSPLLLFLPLLLLLDVP TAAVQASPLQALDFFGNGPP VNYKTGNLYLRGPLKKSNA LVNVTLYYEALCGGCRAFLI RELFPTWLLVMEILNVTLP YGNAQEQNVSGRWEFKCQH EEECKFNKVEACVLDELME LAFLTIVCMEEFEDMERSLP LCLQLYAPGLSPDTIMECAM GDRGMQLMHANAQRD DALQP PHEYVPWVTVNGKPLEDQTQ LLTLVCQLYQGKKPDVCPS
	Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human IFI30 (NP_006323.2).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using IFI30 Rabbit polyclonal antibody (STJ114964) at 1:3000 dilution. _Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. _Lysates/proteins: 25 Mu g per lane. _Blocking buffer: 3% nonfat dry milk in TBST. _Detection: ECL Basic Kit. _Exposure time: 30s.