

Anti-DYNC2H1 antibody (1-250) (STJ115092)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 493kDa Observed Mw: 492kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79659
Gene Symbol	DYNC2H1
UniProt ID	DYHC2_HUMAN
Immunogen Region	1-250
Immunogen Sequence	MANGTADVRKLFITTTQNY FGLMSELWDQPLLCNCLEIN NFLDDGNQMLLRVQRSDAGI SFSNTIEFGDTKDKVLVFFK LRPEVITDENLHDNVLVSSM LESPISSLYQAVRQVFAPML LKDQEWSRNFDPKLNLLSE LEAGLGIVLRSDTNLTCLK FKEDDTRGILTPSDEFQFWI EQAHRGNKQISKERANYFKE LFETIAREFYNLDSLLELV VDLVETTQDVVDDVWRQTE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DYNC2H1 (NP_001073932.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of lysates from rat testis, using DYNC2H1 Rabbit polyclonal antibody (STJ115092) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.