

## Anti-TTC5 antibody (1-300) (STJ115198)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 49kDa Observed Mw: 49kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">91875</a>
Gene Symbol	<a href="#">TTC5</a>
UniProt ID	<a href="#">TTC5_HUMAN</a>
Immunogen Region	1-300
Immunogen Sequence	MMADDEEEVKPILQKLOELV DQLYSFRDCYFETHSVEDAG RKQQDVQKEMEKTLQQMEEV VGSVQGKAQVLMMLTGKALNV TPDYSPKAEELLKAVKLEP ELVEAWNQLGEVYWKKGDDVA AAHTCFSGALHCRNKVSLQ NLSMVLRLRQRTDTEDEHSHH VMDSVROAKLAVQMDVHDGR SWYILGNSYLSLYFSTGQNP KISQQALSAYAQAQEKVDRKA SSNPDLHLNRATLHKYEEES
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TTC5 (NP_612385.2).

### ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates using TTC5 Rabbit polyclonal antibody (STJ115198) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.