

Anti-VGLL3 antibody (71-180 aa) (STJ115228)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 36kDa Observed Mw: 36kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	389136
Gene Symbol	VGLL3
UniProt ID	VGLL3_HUMAN
Immunogen Region	71-180 aa
Immunogen Sequence	EEEEKDQPAEMEYLNSRCVL FTYFQGDIGSVVDEHFSRAL GQAITLHPESAISKSKMGLT PLWRDSSALSSQRNSFPTSF WTSSYQPPAPCLGGVHPDF QVTGPPGTFS
Specificity	A synthetic peptide corresponding to a sequence within amino acids 71-180 of human VGLL3 (NP_057290.2).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates, using VGLL3 rabbit polyclonal antibody (STJ115228) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.