

Anti-LHX8 antibody (1-220) (STJ115416)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 39kDa Observed Mw:
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	431707
Gene Symbol	LHX8
UniProt ID	LHX8_HUMAN
Immunogen Region	1-220
Immunogen Sequence	MQILSRCQGLMSEECGRITTA LAAGRTRKGAGEEGLVSPEG AGDEDCSSSAPLSPSSSPR SMASGSGCPPGKVCVNSCGL EIVDKYLLKVNDLCWHVRCL SCSVCRTSLGRHTSCYIKDK DIFCKLDYFRRYGTRCSRSG RHHSTDWVRRAKGNVYHLA CFACFSCKRQLSTGEEFALV EEKVLCRVHYDCMLDNLKRE VENGNGISVEGALLTEQDVN
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human LHX8 (NP_001001933.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunofluorescence analysis of HeLa cells using LHX8 antibody (STJ115416).