

## Anti-LIG3 antibody (254-460) (STJ115468)

### GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

### PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:1000 IHC-P:1:50-1:200 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 113kDa Observed Mw: 105kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION


Gene ID	<a href="#">3980</a>
Gene Symbol	<a href="#">LIG3</a>
UniProt ID	<a href="#">DNL13_HUMAN</a>
Immunogen Region	254-460
Immunogen Sequence	CDPRHKDCLLREFRKLCAV ADNPSYNTKTIQDFLRKG SAGDGFHGDVYLTVKLLLP VIKTVYNLNDKQIVKLFSAI FNCNPDDMARDLEQGDVSET IRVFFEQSKSFPPAAKSLT IQEVDFLLRLSKLTKEDEQ QQALQDIASRCTANDLKCI RLIKHDLKMNSGAKHVLDAL DPNAYEAFKASRNLDQDVVER VLHNAQE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 254-460 of human DNA Ligase III/LIG3 (NP_039269.2).

### ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.




Western blot analysis of extracts of various cell lines, using (STJ115468) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.




Immunohistochemistry analysis of paraffin-embedded mouse kidney using DNA Ligase III/LIG3 rabbit polyclonal antibody (STJ115468) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.




Immunohistochemistry analysis of paraffin-embedded rat brain using DNA Ligase III/LIG3 rabbit polyclonal antibody (STJ115468) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.




Immunohistochemistry analysis of paraffin-embedded rat liver using DNA Ligase III/LIG3 rabbit polyclonal antibody (STJ115468) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of paraffin-embedded rat ovary using DNA Ligase III/LIG3 rabbit polyclonal antibody (STJ115468) at dilution of 1:50 (40x lens). Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with immunohistochemistry staining protocol.



Immunofluorescence analysis of NIH/3T3 cells using DNA Ligase III/LIG3 rabbit polyclonal antibody (STJ115468) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.



Immunofluorescence analysis of PC-12 cells using DNA Ligase III/LIG3 rabbit polyclonal antibody (STJ115468) at dilution of 1:50 (40x lens). Blue: DAPI for nuclear staining.