

Anti-PSMA4 antibody (1-261) (STJ115496)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat/Monkey

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 29kDa Observed Mw: 29kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5685
Gene Symbol	PSMA4
UniProt ID	PSA4_HUMAN
Immunogen Region	1-261
Immunogen Sequence	MSRRYDSRTTIFSGEGRLYQ VEYAMEAIGHAGTCLGILAN DGVLLAAERRNIHKLLDEVF FSEKIYKLNEDMACSVAGIT SDANVLTNELRLIAQRYLLQ YQEPIPCQLVTALCDIKQA YTQFGGKRPFVSLLYIGWD KHYGFQLYQSDPSGNYGGWK ATCIGNNSAAAVSMLKQDYK EGEMTLKSALALAIKVLNKT MDVSKLSAEKVEIATLTREN GKTIVIRVLKQKEVEQLIKK
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-261 of human PSMA4 (NP_001096137.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of various cell lines, using PSMA4 antibody (STJ115496) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.