

Anti-VEGFD antibody (20-220) (STJ115576)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA: Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 40kDa Observed Mw: 54kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2277
Gene Symbol	VEGFD
UniProt ID	VEGFD_HUMAN
Immunogen Region	20-220
Immunogen Sequence	QGSSNEHGPKVRSSQSTLER SEQQIRAASSLEELLRITHS EDWKLWRCRLRLKSFTSMDS RSASHRSTRFAATFYDIETL KVIDEEWQRTQCSPRETCVE VASELGKSTNTFFKPPCVNV FRCGGCCNEESLICMNTSTS YISKQLFEISVPLTSVPELV PVKVANHTGCKCLPTAPRHP YSIIRRSIQIPEEDRCSHSK K
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-220 of human VEGF-D (NP_004460.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of 293T cells, using VEGFD antibody (STJ115576) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.