

Anti-DPPA4 antibody (15-304) (STJ115677)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 34kDa Observed Mw: 38kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55211
Gene Symbol	DPPA4
UniProt ID	DPPA4_HUMAN
Immunogen Region	15-304
Immunogen Sequence	KGKEWTSTEKSREEDQQASN QPNSIALPGTSAKRTKEKMS IKGSKVLCPPKKKAHEHTDNPR PQKKIPIPPLPSKLPVNL I HRDILRAWCQQLKLSKGGQK LDAYKRLCAFAYPNQKDFPS TAKEAKIRKSLQKKLKVKEG ETSLQSSETHPPEVALPPVG EPPALENSTALLEGVNTVVV TTSAPALLASWARISARAR TPEAVESPQEASGVRWCVVH GKSLPADTDGWVHLQFHAG
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 15-304 of human DPPA4 (NP_060659.3).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates using DPPA4 Rabbit polyclonal antibody (STJ115677) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.