

Anti-NGF antibody (19-241) (STJ115857)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 27kDa Observed Mw: 27kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4803
Gene Symbol	NGF
UniProt ID	NGF_HUMAN
Immunogen Region	19-241
Immunogen Sequence	EPHSESNVPAGHTIPQAHWT KLQHSLDTALRRARSAPAAA IAARVAGQTRNITVDPRLFK KRRLRSPRVLFSTQPPREAA DTQDLDFEVGGAAPFNRTHR SKRSSHPIFHRGEFSVCDS VSVVVGDKTTATDIKIGKEVM VLGEVNNNSVFKQYFFETK CRDPNPVDSGCRGIDSKHWN SYCTTTHTFVKALTMGKQA AWRFIRIDTACVCVLSRKAV RRA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 19-241 of human NGF (NP_002497.2).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of lysates from mouse heart, using NGF Rabbit polyclonal antibody (STJ115857) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS00856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST.