

Anti-EWSR1 antibody (1-300) (STJ115913)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 68kDa Observed Mw:
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2130
Gene Symbol	EWSR1
UniProt ID	EWS_HUMAN
Immunogen Region	1-300
Immunogen Sequence	MASTDYSTYSQAAAQQGYSA YTAQPTQGYAQTQAYGQQS YGTYGQPTDVSYTQAQTAT YGQTAYATSYGQPPTGYTTP TAPQAYSQPVQGYGTGAYDT TTATVTTTQASYAAQSAYGT QPAYPAYGQQPAATAPTRPQ DGNKPTETSQPSSTGGYNQ PSLGYGQSNYSYPQVPGSYP MQPVTAPPSYPPTSYSSTQP TSYDQSSYSQQNTYGQPSSY GQQSSYGQQSSYGQQPPTS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human EWSR1 (NP_001156757.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunofluorescence analysis of U2OS cells using EWSR1 antibody (STJ115913).