

Anti-SETD5 antibody (520-740) (STJ116075)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 158kDa Observed Mw:
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55209
Gene Symbol	SETD5
UniProt ID	SETD5_HUMAN
Immunogen Region	520-740
Immunogen Sequence	HAFENLEKRKKRRDQPLEQS NSDVEITTTTSETPVGEETK TEAPESEVSNVSNSVTIPST PQSVGVNTRRSSQAGDIAAE KLVKPPPAKPSRPRPKSRI SRYRTSSAQLRKQKQANAQ QAELSQAALIEGGNSLVTP TEAGSLDSSGENRPLTGSDP TVVSITGSHVNRAASKYPKT KKYLVTEWLNDAEKQCEPV ECPLRITTDPTVLATTLNML P
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 520-740 of human SETD5 (NP_001073986.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunofluorescence analysis of MCF-7 cells using SETD5 antibody (STJ116075). Blue: DAPI for nuclear staining.