

Anti-GSTO2 antibody (124-243) (STJ116091)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IF/ICC:1:50-1:100 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 28kDa Observed Mw:
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	119391
Gene Symbol	GSTO2
UniProt ID	GSTO2_HUMAN
Immunogen Region	124-243
Immunogen Sequence	PHLTKECLVALRCGRECTNL KAALRQEFSNLEEILEYQNT TFFGGTCISMIDYLLWPWFE RLDYYGILDCVSHTPALRLW ISAMKWDPTVCALLMDKSIF QGFLNLYFQNNPNAFDGLC
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 124-243 of human GSTO2 (NP_899062.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Immunofluorescence analysis of U2OS cells using GSTO2 antibody (STJ116091). Blue: DAPI for nuclear staining.