

Anti-GGT5 antibody (388-587) (STJ116586)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 62kDa Observed Mw: 62kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2687
Gene Symbol	GGT5
UniProt ID	GGT5_HUMAN
Immunogen Region	388-587
Immunogen Sequence	TSHVSVLGEDGSAVAATSTI NTPFGAMVYSPRTGIILNNE LLDLCERCPRGSGTTPSPAV SGDRVGGAPGRCWPPVPGER SPSSMVPSILINKAQGSKLV IGGAGGELIISAVAQAIMSK LWLGFDLRAAIAAPILHVNS KGCVEYEPNFSQEVQRGLQD RGQNQTQRPFFLNVVQAVSQ EGACVYAVSDLRKSGEAGY
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 388-587 of human GGT5 (NP_001093251.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of lysates from HT-29 cells, using GGT5 Rabbit polyclonal antibody (STJ116586) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.