

Anti-MIOX antibody (1-285) (STJ116619)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/ICC:1:50-1:200 ELISA:Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 33kDa Observed Mw: 33kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	55586
Gene Symbol	MIOX
UniProt ID	MIOX_HUMAN
Immunogen Region	1-285
Immunogen Sequence	MKVTVGPDPSLVYRPDVPDPE VAKDKASFRNYTSGPLLDRV FTTYKLMHHTQTVDFVRSKH AQFGGFSYKKMTVMEAVDLL DGLVDESDDPDVDFPNSFHAF QTAEGIRKAHPDKDWFHLVG LLHDLGKVLALFGEPQWAVV GDTFPVGCRRPQASVVFCDST FQDNPDLQDPRYSTELGMYQ PHCGLDRVLMSWGHDYMYQ VMKFNKFSLPPEAFYMIRFH SFYPWHTGRDYQQLCSSQQD
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-285 of human MIOX (NP_060054.4).

ADDITIONAL INFORMATION


Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of paraffin-embedded Mouse kidney using MIOX Rabbit polyclonal antibody (STJ116619) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of paraffin-embedded Rat kidney using MIOX Rabbit polyclonal antibody (STJ116619) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using MIOX Rabbit polyclonal antibody (STJ116619) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.