

Anti-B3GNT5 antibody (36-378) (STJ116638)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 44kDa Observed Mw: 44kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84002
Gene Symbol	B3GNT5
UniProt ID	B3GN5_HUMAN
Immunogen Region	36-378
Immunogen Sequence	DNHIVSHMKSYRYRLINSY DFVNDTSLKHTSAGPRYQY LINHKEKCAQDVLVLLLFVK TAPENYDRRSGIRRTWGNEN YVRSQNLNANIKTLFALGTPN PLEGEELQRKLAWEDQRYND IIQQDFVDSFYNLTLKLLMQ FSWANTYCPHAKFLMTADD IFIHMPNLIEYLOSLEQIGV QDFWIGRVHRGAPPIRDKSS KYVVSYEMYQWPAYPDYTAG AAYVISGDVAAKVYEASQT
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 36-378 of human B3GNT5 (NP_114436.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates using B3GNT5 Rabbit polyclonal antibody (STJ116638) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.