

Anti-C7orf50 antibody (1-194) (STJ116694)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 22kDa Observed Mw: Refer to Figures
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84310
Gene Symbol	CHLSN
UniProt ID	CHOLN_HUMAN
Immunogen Region	1-194
Immunogen Sequence	MAKQKRKVPVEVTEKKNKKLK KASAEGPLLGPEAAPSGEA GSKGEAVLRPGLDAEPELSP EEQVLERKLLKKERKKEERQ RLREAGLVAQHPPARRSGAE LALDYLCRWAKHKNWRFQK TRQTWLLHMYDSKVPDEH FSTLLAYLEGLQGRARELTV QKAEALMRELDEEGSDPPLP GRAQRIQVLQLLS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human C7orf50 (NP_115726.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of various cell lines, using C7orf50 antibody (STJ116694) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.