

Anti-ZAK antibody (100-312) (STJ116724)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 35kDa/51kDa/91kDa Observed Mw: 105kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	51776
Gene Symbol	MAP3K20
UniProt ID	M3K20_HUMAN
Immunogen Region	100-312
Immunogen Sequence	EEMDMDHIMTWDVAKGMH YLHMEAPVKVIHRDLKSRNV VIAADGVLKICDFGASRFHN HTTHMSLVGTFPWMAPEVIQ SLPVSETCDTYSYGVVLWEM LTREVPFKGLEGLQVAVLVV EKNERLTIPSSCPRSFAELL HQCWEADAKKRPSFKQIISI LESMSNDTSLPDKCNSFLHN KAEWRCEIEATLERLKKLER DLSFKQEQLKERE
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 100-312 of human ZAK (NP_598407.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of L929 cells using ZAK Rabbit polyclonal antibody (STJ116724) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using ZAK Rabbit polyclonal antibody (STJ116724) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.