

Anti-JUNB antibody (1-200) (STJ116777)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:100-1:500 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 36kDa Observed Mw: 42kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3726
Gene Symbol	JUNB
UniProt ID	JUNB_HUMAN
Immunogen Region	1-200
Immunogen Sequence	MCTKMEQPFYHDDSYTATGY GRAPGGLSLHDYKLLKPSLA VNLADPYRSLKAPGARGPGP EGGGGGSYFSGQSDTGASL KLASSELERLIVPNSNGVIT TTPPPGQYFYPRGGGSGGG AGGAGGGVTEEQEGFADGFV KALDDLHKMNHVTPPNVSLG ATGGPPAGPGGVYAGPEPPP VYTNLSSYSPASASSGGAGA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human JunB (NP_002220.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of various lysates using JunB Rabbit polyclonal antibody (STJ116777) at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.