

Anti-SNRNP70 antibody (1-180) (STJ116986)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 52kDa Observed Mw: 63kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6625
Gene Symbol	SNRNP70
UniProt ID	RU17_HUMAN
Immunogen Region	1-180
Immunogen Sequence	MTQFLPPNLLALFAPRDPIP YLPPLEKLPHEKHHNQPYCG IAPYIREFEDPRDAPPPTRA ETREERMERKRREKIERRQQ EVETELKMWDPHNDPNAQGD AFKTLFVARVNYDTTESKLR REFEVYGPRIKRIHMYYSKRS GKPRGYAFIEYEHERDMHSA YKHADGKKIDGRRVLVDVER
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human SNRNP70 (NP_003080.2).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of HeLa, using SNRNP70 rabbit polyclonal antibody (STJ116986) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% non-fat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.