

Anti-SLC16A4 antibody (388-487) (STJ117018)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 54kDa Observed Mw: 60kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9122
Gene Symbol	SLC16A4
UniProt ID	MOT5_HUMAN
Immunogen Region	388-487
Immunogen Sequence	TFPLLMTYTICFAIFAGGYL ALILPVLVDLCRNSTVNRFL GLASFFAGMAVLSGPPIAGW LYDYTQTYNGSFYFSGICYL LSSVSFFFVPLAERWKNSLT
Specificity	A synthetic peptide corresponding to a sequence within amino acids 388-487 of human SLC16A4 (NP_004687.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of extracts of SW480 cells, using SLC16A4 Rabbit polyclonal antibody (STJ117018) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.