

Anti-FBXO42 antibody (468-717) (STJ117098)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 78kDa Observed Mw: 78kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	54455
Gene Symbol	FBXO42
UniProt ID	FBX42_HUMAN
Immunogen Region	468-717
Immunogen Sequence	PSAPEGYDLKIGLSLAPRRG SLPDQKDLRLGSDLNWDLK PASSSNPMDGMDNRTVGGSM RHPPEQTNGVHTPPHVASAL AGAVSPGALRRSLEAIKAMS SKGPSASAALSPPLGSSPGS PGSQSLSSGETVPIPRPGPA QGDGHSLLPIARRLGHHPPO SLNVGKPLYQSMNCKPMQMY VLDIKDTKEKGRVKWKVFNS SSVVGPPETSLHTVVQGRGE LIIFGGLMDKKQNVKYYPK
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 468-717 of human FBXO42 (NP_061867.1).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of mouse lung, using FBXO42 antibody (STJ117098) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.