

Anti-B3GNT6 antibody (30-200) (STJ117174)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 43kDa Observed Mw: 43kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	192134
Gene Symbol	B3GNT6
UniProt ID	B3GN6_HUMAN
Immunogen Region	30-200
Immunogen Sequence	FLQAPRSPREERSPQEETPE GPTDAPAADEPPSELVPGPP CVANASANATADFEQLPARI QDFLRVYRHCRRHFLLWDAPA KCAGGRGVFLLLVKSAPEH YERRELIRRTWQERSYGGGR PVRRLFLLTGTPGPEDEARAE RLAELVALEAREHGDVLQWA FADTFNLNLTk
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 30-200 of human B3GNT6 (NP_619651.3).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of lysates from rat stomach, using B3GNT6 Rabbit polyclonal antibody (STJ117174) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 60s.