

Anti-DAPK3 antibody (1-100) (STJ117241)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IHC-P/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:100 IF/CC:1:50-1:100 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 53kDa Observed Mw: 55kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	1613
Gene Symbol	DAPK3
UniProt ID	DAPK3_HUMAN
Immunogen Region	1-100
Immunogen Sequence	MSTFRQEDVEDHYEMGEELG SGQFAIVRKCRQKGTGKEYA AKFIKKRRLSSRRGVSREE IEREVNLRIRHPNIITLH DIFENKTDVVLILELVSGGE
Specificity	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DAPK3 (NP_001339.1).

ADDITIONAL INFORMATION


Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**




Immunofluorescence analysis of U-2 OS cells using DAPK3 Rabbit polyclonal antibody (STJ117241) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.




Immunofluorescence analysis of C6 cells using DAPK3 Rabbit polyclonal antibody (STJ117241) at dilution of 1:100 (40x lens). Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.




Immunohistochemistry analysis of DAPK3 in paraffin-embedded mouse kidney using DAPK3 Rabbit polyclonal antibody (STJ117241) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.




Immunohistochemistry analysis of DAPK3 in paraffin-embedded human placenta using DAPK3 Rabbit polyclonal antibody (STJ117241) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of DAPK3 in paraffin-embedded human lung cancer using DAPK3 Rabbit polyclonal antibody (STJ117241) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.



Immunohistochemistry analysis of DAPK3 in paraffin-embedded rat kidney using DAPK3 Rabbit polyclonal antibody (STJ117241) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7. 2 before commencing with immunohistochemistry staining protocol.



Western blot analysis of various lysates using DAPK3 Rabbit polyclonal antibody (STJ117241) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.