

Anti-DCAF11 antibody (1-170) (STJ117714)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:1000-1:5000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 62kDa Observed Mw: 61kDa
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80344
Gene Symbol	DCAF11
UniProt ID	DCA11_HUMAN
Immunogen Region	1-170
Immunogen Sequence	MGSRNSSSAGSGSGDPSEGL PRRGAGLRSEEEEEDEDV DLAQVLAYLLRRGQVRLVQG GGAANLQFIQALLDSEEND RAWDGRLGDRYNPPVDATPD TRELEFNEIKTQVELATGQL GLRRAAQKHSFPRMLHQRER GLCHRGSFSLGGEQSRVISHF LPNDLGFTDS
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human DCAF11 (NP_079506.3).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of Mouse liver, using DCAF11 antibody (STJ117714) at 1:2000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 20s.