

Anti-ZNF598 antibody (120-420) (STJ117741)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:200-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 99kDa Observed Mw: 115kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	90850
Gene Symbol	ZNF598
UniProt ID	ZNF598_HUMAN
Immunogen Region	120-420
Immunogen Sequence	PRCPPELPPFSLFGDLEQHMRRQHELFCRCRLCLQLHLQIFTYERKWYSRKDLARHRMQGDPD DTSHRGHPLCKFCDERLYLDNDELLKHLRRDHYFCHFCDSDGAQDYYSYAYLREHFREKH FLCEEGRCSSTEQFTHAFRTEIDLKAHRTACHSRRAEARQNRHIDLQFSYAPRHSRRNEG VVGGEDYEEVDYRSRQGRVAVRAGTRGAQQSRRGSRWRYKREEDREVAAVRASVAAQQQ
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 120-420 of human ZNF598 (NP_835461.2).

ADDITIONAL INFORMATION

Note **STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.**



Western blot analysis of various lysates using ZNF598 Rabbit polyclonal antibody (STJ117741) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15S.