

Anti-TNFRSF9 antibody (164-255) (STJ117796)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/IHC-P/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:500-1:2000 IHC-P:1:50-1:200 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH 7.3.
	Isotype	IgG
	Molecular Weight	Protein Mw: 28kDa Observed Mw: 27kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	3604
	Gene Symbol	TNFRSF9
	UniProt ID	TNR9_HUMAN
	Immunogen Region	164-255
	Immunogen Sequence	SMCEQDCKQGQELTKKGCKD CCFGTFNDQKRGICRPWTNC SLDGKSVLVNGTKERDVVCG PSPADLSPGASSVTPAPAR EPGHSPQIISFFLALTSTAL LFLFFLTLRFSVVKRGRKK LLYIFKQPFMRPVQTTQEED GCSCRFPEEEEGGCEL
	Specificity	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human 4-1BB/CD137 (NP_001552.2).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
   		