

Anti-IL10RB antibody (20-220) (STJ118126)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 IF/CC:1:50-1:200 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.05% Proclin300, 50% Glycerol, pH7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 37kDa Observed Mw: 36kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	3588
Gene Symbol	IL10RB
UniProt ID	I10R2_HUMAN
Immunogen Region	20-220
Immunogen Sequence	MVPPPENVRMNSVNFKNILQ WESPAFAKGNLTFTAQYLSY RIFQDKCMNTTLTECDFSSL SKYGDHTLRVRAEFADEHSD WVNITFCPVDDTIIGPPGMQ VEVLADSLHMRFLAPKIENE YETWTMKNVYNSWTYNVQYW KNGTDEKFQITPQYDFEVLRL NLEPWTTYCVQVRGFLPDRN KAGEWSEPVCEQTTHDETVP S
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 20-220 of human IL10RB (NP_000619.3).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of L929 cells using IL10RB Rabbit polyclonal antibody (STJ118126) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) at 1:500 dilution. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using IL10RB Rabbit polyclonal antibody (STJ118126) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.