

Anti-RBP7 antibody (1-134) (STJ118384)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/IF/ICC/ELISA
Host / Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA: Recommended starting concentration is 1 μ g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 16kDa Observed Mw: 15kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	116362
Gene Symbol	RBP7
UniProt ID	RET7_HUMAN
Immunogen Region	1-134
Immunogen Sequence	MPADLSGTWTLSSDNFEGY MLALGIDFATRKIALLKPQ KVIEQNGDSFTIHTNSSLRN YFVKFKVGEEFDEDNRGLDN RKCKSLVIWDNDRLTCIQKG EKKNRGWTHWIEGDKLHLEM FCEGQVCKQTFQRA
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-134 of human RBP7 (NP_443192.1).

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using RBP7 Rabbit polyclonal antibody (STJ118384) at 1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 μ g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.