


Anti-INTS7 antibody (702-948) (STJ118553)

GENERAL INFORMATION		
	Product Type	Primary antibodies
	Applications	WB/ELISA
	Host / Source	Rabbit
	Reactivity	Human/Mouse
PRODUCT PROPERTIES		
	Clonality	Polyclonal
	Concentration	Lot specific
	Conjugation	Unconjugated
	Purification	Affinity purification
	Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
	Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
	Isotype	IgG
	Molecular Weight	Protein Mw: 107kDa Observed Mw: 107kDa
	Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Gene ID	25896
	Gene Symbol	INTS7
	UniProt ID	INT7_HUMAN
	Immunogen Region	702-948
	Immunogen Sequence	NVELQQQSCLLISHAIEALI LDPESASFQEYGSTGTAHAD SEYERRMMSVYNHVL EEHTA CLCNAILALKVPLSFQRYF FQKLQSTSIKLALSPSPRNP AEPIAVQNNQQLALKVEGVV QHGSKPGLFRKIQSVCLNVS STLQSKSGQDYKIPIDNMTN EMEQRVEPHNDYFSTQFLLN FAILGTHNITVESSVKDANG IVWKTGPRTTIFVKSLEDPY SQQIRLQQQQAQQPLQQQQ
	Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 702-948 of human INTS7 (NP_001186740.1).
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Western blot analysis of various lysates using INTS7 Rabbit polyclonal antibody (STJ118553) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 Mu g per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.