

## Anti-SGOL1 antibody (10-100) (STJ118614)

### GENERAL INFORMATION

|               |                    |
|---------------|--------------------|
| Product Type  | Primary antibodies |
| Applications  | WB/IF/ICC/ELISA    |
| Host / Source | Rabbit             |
| Reactivity    | Human/Mouse/Rat    |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                                          |
|---------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Clonality           | Polyclonal                                                                                                                                                               |
| Concentration       | Lot specific                                                                                                                                                             |
| Conjugation         | Unconjugated                                                                                                                                                             |
| Purification        | Affinity purification                                                                                                                                                    |
| Dilution Range      | WB:1:500-1:2000<br>IF/CC:1:50-1:200<br>ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements. |
| Formulation         | PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.                                                                                                                         |
| Isotype             | IgG                                                                                                                                                                      |
| Molecular Weight    | Protein Mw: 64kDa<br>Observed Mw: 75kDa                                                                                                                                  |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                                                                           |

### TARGET INFORMATION


|                    |                                                                                                                   |
|--------------------|-------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">151648</a>                                                                                            |
| Gene Symbol        | <a href="#">SGO1</a>                                                                                              |
| UniProt ID         | <a href="#">SGO1_HUMAN</a>                                                                                        |
| Immunogen Region   | 10-100                                                                                                            |
| Immunogen Sequence | SFQDSLEDIKKRMKEKRNKN LAEIGKRRSFIAAPCQIITN TSTLLKNYQDNNKMLVLALE<br>NEKSKVKEAQDIIQLRKEC YYLTCQLYALK                 |
| Specificity        | Recombinant fusion protein containing a sequence corresponding to amino acids 10-100 of human SGO1 (NP_612493.1). |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Immunofluorescence analysis of U2OS cells using SGOL1 Rabbit polyclonal antibody (STJ118614) at dilution of 1:100. Blue: DAPI for nuclear staining.



Western blot analysis of various lysates using SGOL1 Rabbit polyclonal antibody (STJ118614) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.