

Anti-ENDOD1 antibody (22-330) (STJ118929)

GENERAL INFORMATION	
Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Human/Mouse/Rat
PRODUCT PROPERTIES	
Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 Mu g/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.01% Thimerosal, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 55kDa Observed Mw: Refer to figures
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	23052
Gene Symbol	ENDOD1
UniProt ID	ENDD1_HUMAN
Immunogen Region	22-330
Immunogen Sequence	RLVGEEEGAGFGECDKFFYAG TPPAGLAADSHVKICQRAEG AERFATLYSTRDRIPVYSAF RAPRPAPGGAEQRWLVEPQI DDPNSNLEEAINAEAITSV NSLGSKQALNTDYLDSDYQR GQLYPFSLSSDVQVATFTLT NSAPMTQSFQERWYVNLHSL MDRALTPQCGSGEDLYILTG TVPSDYRVKDKVAVPEFVWL AACCAVPGGGWAMGFVKHTR DSDIIEDVMVKDLQKLLPF
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 22-330 of human ENDOD1 (NP_055851.1).
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.