

Anti-NLRC5 antibody (1-240) (STJ119122)

GENERAL INFORMATION

Product Type	Primary antibodies
Applications	WB/ELISA
Host / Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Concentration	Lot specific
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB:1:500-1:2000 ELISA:Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.
Formulation	PBS with 0.09% Sodium Azide, 50% Glycerol, pH 7.3.
Isotype	IgG
Molecular Weight	Protein Mw: 205kDa Observed Mw: 202kDa
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84166
Gene Symbol	NLRC5
UniProt ID	NLRC5_HUMAN
Immunogen Region	1-240
Immunogen Sequence	MDPVGLQLGNKNLWSCLVRL LTKDPEWLNAMKFFLPNTD LDSRNETLDPEQRVILQLNK LHVQGS DWTWQSFHCVCMQL EVPLDLEVL LSTFGYDDGF TSQLGAEGKSQPESQLHHGL KRPHQSCGSSPRRKQCKKQQ LELAKKYLQLLRTSAQQRYR SQIPGSGQPHAFHQVYVPI LRRATASLDTPEGAIMGDVK VEDGADVSISDLFNTRVNKG PRVTVLLGKAGMGKTTLAH
Specificity	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human NLRC5 (NP_115582.4).

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



Western blot analysis of extracts of mouse lung, using NLRC5 antibody (STJ119122) at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (STJS000856) at 1:10000 dilution. Lysates/proteins: 25 µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5min.